


ANALISIS GRANULOMETRICOS

MUESTRA: D744

UTM X: 645509


UTM Y: 3244525



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.36 (100 %C)	0.4%	99.6%
n°10	2 mm	1.46 (100 %C)	1.8%	98.2%
n°18	1 mm	12.15 (100 %C)	14.0%	86.0%
n°25	0.71 mm	16.28 (100 %C)	30.3%	69.7%
n°35	0.5 mm	19.81 (0 %C)	50.1%	49.9%
n°60	0.25 mm	34.63 (0 %C)	84.7%	15.3%
n°80	0.18 mm	8.38 (0 %C)	93.1%	6.9%
n°120	0.125 mm	5.47 (0 %C)	98.5%	1.5%
n°230	0.063 mm	1.09 (0 %C)	99.6%	0.4%
Finos	<0,063 mm	0.36 (0 %C)	100.0%	0.0%

% FINOS	0.52%
% CASCAJO	30.26%
D16	0.592
D50	0.391
D84	0.215
Media	0.399
σ	0.188
M. distribuc.	0.404

MODA ARENAS MEDIAS



	<i>D50</i>	
<i>Cantos y gravas</i>	>2	mm
<i>Arena muy gruesa</i>	1-2	mm
<i>Arena gruesa</i>	0,5-1	mm
<i>Arena media</i>	0,25-0,5	mm
<i>Arena fina</i>	0,125-0,25	mm
<i>Arena muy fina</i>	0,062-0,125	mm
<i>Fango</i>	<0,062	mm