


ANALISIS GRANULOMETRICOS

MUESTRA: D761

UTM X: 652192


UTM Y: 3222949



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.38 (100 %C)	0.4%	99.6%
n°10	2 mm	4.53 (100 %C)	4.9%	95.1%
n°18	1 mm	17.09 (100 %C)	22.0%	78.0%
n°25	0.71 mm	12.46 (100 %C)	34.5%	65.5%
n°35	0.5 mm	9.35 (0 %C)	43.8%	56.2%
n°60	0.25 mm	22.47 (0 %C)	66.3%	33.7%
n°80	0.18 mm	8.59 (0 %C)	74.9%	25.1%
n°120	0.125 mm	14.16 (0 %C)	89.0%	11.0%
n°230	0.063 mm	9.73 (0 %C)	98.8%	1.2%
Finos	<0,063 mm	1.23 (0 %C)	100.0%	0.0%

% FINOS	1.87%
% CASCAJO	34.47%
D16	0.487
D50	0.242
D84	0.122
Media	0.284
σ	0.183
M. distribuc.	0.305

MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm