


ANALISIS GRANULOMETRICOS

MUESTRA: IP111

UTM X: 633145


UTM Y: 3200658



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.12 (50 %C)	0.1%	99.9%
n°18	1 mm	1.17 (50 %C)	1.3%	98.7%
n°25	0.71 mm	2.21 (25 %C)	3.5%	96.5%
n°35	0.5 mm	5.29 (0 %C)	8.8%	91.2%
n°60	0.25 mm	36.68 (0 %C)	45.5%	54.5%
n°80	0.18 mm	32.60 (0 %C)	78.1%	21.9%
n°120	0.125 mm	20.23 (0 %C)	98.3%	1.7%
n°230	0.063 mm	1.71 (0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.00 (0 %C)	100.0%	0.0%

% FINOS	0.00%
% CASCAJO	1.20%
D16	0.444
D50	0.239
D84	0.163
Media	0.282
σ	0.140
M. distribuc.	0.304

MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm