


ANALISIS GRANULOMETRICOS

MUESTRA: IP093

UTM X: 613949


UTM Y: 3204057



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.09 (0 %C)	0.1%	99.9%
n°18	1 mm	0.26 (5 %C)	0.3%	99.7%
n°25	0.71 mm	1.90 (5 %C)	2.2%	97.8%
n°35	0.5 mm	8.04 (0 %C)	10.3%	89.7%
n°60	0.25 mm	84.62 (0 %C)	94.9%	5.1%
n°80	0.18 mm	4.24 (0 %C)	99.1%	0.9%
n°120	0.125 mm	0.69 (0 %C)	99.8%	0.2%
n°230	0.063 mm	0.09 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.09 (0 %C)	100.0%	0.0%

% FINOS	0.09%
% CASCAJO	0.11%
D16	0.483
D50	0.383
D84	0.282
Media	0.383
σ	0.100
M. distribuc.	0.383

MODA ARENAS MEDIAS



	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm