


ANALISIS GRANULOMETRICOS

MUESTRA: IP103

UTM X: 634147


UTM Y: 3200977



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.00 (100 %C)	0.0%	100.0%
n°25	0.71 mm	0.04 (50 %C)	0.0%	100.0%
n°35	0.5 mm	0.13 (50 %C)	0.2%	99.8%
n°60	0.25 mm	7.04 (0 %C)	7.2%	92.8%
n°80	0.18 mm	27.22 (0 %C)	34.4%	65.6%
n°120	0.125 mm	59.10 (0 %C)	93.5%	6.5%
n°230	0.063 mm	6.46 (0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.00 (0 %C)	100.0%	0.0%

% FINOS	0.00%
% CASCAJO	0.09%
D16	0.227
D50	0.165
D84	0.134
Media	0.176
σ	0.047
M. distribuc.	0.181

MODA ARENAS FINAS



	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm