


## ANALISIS GRANULOMETRICOS

MUESTRA: IP48

UTM X: 630956


UTM Y: 3221367



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 ( 0 %C)	0.0%	100.0%
n°10	2 mm	0.13 ( 25 %C)	0.1%	99.9%
n°18	1 mm	0.44 ( 75 %C)	0.6%	99.4%
n°25	0.71 mm	0.50 ( 90 %C)	1.1%	98.9%
n°35	0.5 mm	4.20 ( 0 %C)	5.3%	94.7%
n°60	0.25 mm	66.52 ( 0 %C)	71.8%	28.2%
n°80	0.18 mm	19.12 ( 0 %C)	90.9%	9.1%
n°120	0.125 mm	8.03 ( 0 %C)	98.9%	1.1%
n°230	0.063 mm	0.82 ( 0 %C)	99.7%	0.3%
Finos	<0,063 mm	0.25 ( 0 %C)	100.0%	0.0%

% FINOS	0.25%
% CASCAJO	0.81%
<b>D16</b>	0.457
<b>D50</b>	0.330
<b>D84</b>	0.205
<b>Media</b>	0.331
<b>σ</b>	0.126
<b>M. distribuc.</b>	0.331

### MODA ARENAS MEDIAS



	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm