


ANALISIS GRANULOMETRICOS

MUESTRA: IP46

UTM X: 630693


UTM Y: 3221395



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.18 (50 %C)	0.2%	99.8%
n°25	0.71 mm	0.79 (75 %C)	1.0%	99.0%
n°35	0.5 mm	7.63 (0 %C)	8.6%	91.4%
n°60	0.25 mm	61.49 (0 %C)	70.1%	29.9%
n°80	0.18 mm	19.30 (0 %C)	89.4%	10.6%
n°120	0.125 mm	9.56 (0 %C)	98.9%	1.1%
n°230	0.063 mm	0.96 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.09 (0 %C)	100.0%	0.0%

% FINOS	0.09%
% CASCAJO	0.68%
D16	0.468
D50	0.330
D84	0.199
Media	0.332
σ	0.134
M. distribuc.	0.333

MODA ARENAS MEDIAS



	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm