


## ANALISIS GRANULOMETRICOS

MUESTRA: D753

UTM X: 654978


UTM Y: 3226172



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.10 ( 100 %C)	0.1%	99.9%
n°10	2 mm	0.96 ( 100 %C)	1.1%	98.9%
n°18	1 mm	10.30 ( 100 %C)	11.4%	88.6%
n°25	0.71 mm	11.07 ( 100 %C)	22.4%	77.6%
n°35	0.5 mm	12.42 ( 0 %C)	34.8%	65.2%
n°60	0.25 mm	49.18 ( 0 %C)	84.0%	16.0%
n°80	0.18 mm	11.36 ( 0 %C)	95.4%	4.6%
n°120	0.125 mm	4.04 ( 0 %C)	99.4%	0.6%
n°230	0.063 mm	0.38 ( 0 %C)	99.8%	0.2%
Finos	<0,063 mm	0.19 ( 0 %C)	100.0%	0.0%

% FINOS	0.25%
% CASCAJO	22.43%
<b>D16</b>	0.500
<b>D50</b>	0.366
<b>D84</b>	0.228
<b>Media</b>	0.365
<b>σ</b>	0.136
<b>M. distribuc.</b>	0.364

## MODA ARENAS MEDIAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm