


# ANALISIS GRANULOMETRICOS

MUESTRA: IP086


UTM X: 644597

UTM Y: 3206841



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 ( 0 %C)	0.0%	100.0%
n°10	2 mm	0.23 ( 25 %C)	0.2%	99.8%
n°18	1 mm	0.77 ( 50 %C)	1.0%	99.0%
n°25	0.71 mm	9.04 ( 80 %C)	10.0%	90.0%
n°35	0.5 mm	45.94 ( 0 %C)	56.0%	44.0%
n°60	0.25 mm	40.20 ( 0 %C)	96.2%	3.8%
n°80	0.18 mm	3.06 ( 0 %C)	99.2%	0.8%
n°120	0.125 mm	0.61 ( 0 %C)	99.8%	0.2%
n°230	0.063 mm	0.08 ( 0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.08 ( 0 %C)	100.0%	0.0%

% FINOS	0.08%
% CASCAJO	7.67%
D16	0.653
D50	0.510
D84	0.318
Media	0.494
σ	0.168
M. distribuc.	0.486



## MODA ARENAS GRUESAS

	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm