


ANALISIS GRANULOMETRICOS

MUESTRA: IP139


UTM X: 619275

UTM Y: 3215515



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	8.85 (0 %C)	8.8%	91.2%
n°10	2 mm	28.78 (0 %C)	37.6%	62.4%
n°18	1 mm	49.22 (0 %C)	86.8%	13.2%
n°25	0.71 mm	11.90 (0 %C)	98.7%	1.3%
n°35	0.5 mm	1.12 (0 %C)	99.9%	0.1%
n°60	0.25 mm	0.14 (0 %C)	100.0%	0.0%
n°80	0.18 mm	0.00 (0 %C)	100.0%	0.0%
n°120	0.125 mm	0.00 (0 %C)	100.0%	0.0%
n°230	0.063 mm	0.00 (0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.00 (0 %C)	100.0%	0.0%

% FINOS	0.00%
% CASCAJO	0.00%
D16	4.067
D50	1.749
D84	1.058
Media	2.291
σ	1.504
M. distribuc.	2.562



MODA ARENAS MUY GRUESAS

	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm