


ANALISIS GRANULOMETRICOS

MUESTRA: IP81

UTM X: 614142


UTM Y: 3193338



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.07 (0 %C)	0.1%	99.9%
n°25	0.71 mm	0.15 (100 %C)	0.2%	99.8%
n°35	0.5 mm	0.96 (90 %C)	1.2%	98.8%
n°60	0.25 mm	55.96 (0 %C)	57.1%	42.9%
n°80	0.18 mm	32.06 (0 %C)	89.2%	10.8%
n°120	0.125 mm	10.37 (0 %C)	99.6%	0.4%
n°230	0.063 mm	0.37 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.07 (0 %C)	100.0%	0.0%

% FINOS	0.07%
% CASCAJO	1.01%
D16	0.430
D50	0.280
D84	0.191
Media	0.300
σ	0.120
M. distribuc.	0.310

MODA ARENAS MEDIAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm