


ANALISIS GRANULOMETRICOS

MUESTRA: **IP1**

UTM X: **642209**


UTM Y: **3233327**



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.07 (100 %C)	0.1%	99.9%
n°18	1 mm	0.13 (50 %C)	0.2%	99.8%
n°25	0.71 mm	0.86 (50 %C)	1.1%	98.9%
n°35	0.5 mm	4.69 (0 %C)	5.7%	94.3%
n°60	0.25 mm	76.11 (0 %C)	81.8%	18.2%
n°80	0.18 mm	14.52 (0 %C)	96.4%	3.6%
n°120	0.125 mm	3.10 (0 %C)	99.5%	0.5%
n°230	0.063 mm	0.46 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.07 (0 %C)	100.0%	0.0%

% FINOS	0.07%
% CASCAJO	0.56%
D16	0.465
D50	0.354
D84	0.239
Media	0.353
σ	0.113
M. distribuc.	0.352

MODA ARENAS MEDIAS



	<i>D50</i>	
<i>Cantos y gravas</i>	>2	mm
<i>Arena muy gruesa</i>	1-2	mm
<i>Arena gruesa</i>	0,5-1	mm
<i>Arena media</i>	0,25-0,5	mm
<i>Arena fina</i>	0,125-0,25	mm
<i>Arena muy fina</i>	0,062-0,125	mm
<i>Fango</i>	<0,062	mm