


ANALISIS GRANULOMETRICOS

MUESTRA: IP090


UTM X: 613920

UTM Y: 3201457



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	15.21 (0 %C)	15.2%	84.8%
n°10	2 mm	84.62 (0 %C)	99.8%	0.2%
n°18	1 mm	0.00 (0 %C)	99.8%	0.2%
n°25	0.71 mm	0.00 (0 %C)	99.8%	0.2%
n°35	0.5 mm	0.00 (0 %C)	99.8%	0.2%
n°60	0.25 mm	0.04 (0 %C)	99.9%	0.1%
n°80	0.18 mm	0.04 (0 %C)	99.9%	0.1%
n°120	0.125 mm	0.04 (0 %C)	100.0%	0.0%
n°230	0.063 mm	0.04 (0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.00 (0 %C)	100.0%	0.0%

% FINOS	0.00%
% CASCAJO	0.00%
D16	4.724
D50	3.620
D84	2.515
Media	3.620
σ	1.105
M. distribuc.	3.620



MODA CANTOS Y GRAVAS

	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm