


ANALISIS GRANULOMETRICOS

MUESTRA: IP108

UTM X: 633429


UTM Y: 3200836



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.08 (50 %C)	0.1%	99.9%
n°25	0.71 mm	0.24 (50 %C)	0.3%	99.7%
n°35	0.5 mm	0.87 (50 %C)	1.2%	98.8%
n°60	0.25 mm	18.44 (0 %C)	19.6%	80.4%
n°80	0.18 mm	37.40 (0 %C)	57.0%	43.0%
n°120	0.125 mm	37.12 (0 %C)	94.2%	5.8%
n°230	0.063 mm	5.76 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.08 (0 %C)	100.0%	0.0%

% FINOS	0.08%
% CASCAJO	0.60%
D16	0.292
D50	0.193
D84	0.140
Media	0.208
σ	0.076
M. distribuc.	0.216

MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm