


## ANALISIS GRANULOMETRICOS

MUESTRA: IP101

UTM X: 634265


UTM Y: 3201067



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 ( 0 %C)	0.0%	100.0%
n°10	2 mm	0.00 ( 0 %C)	0.0%	100.0%
n°18	1 mm	0.00 ( 100 %C)	0.0%	100.0%
n°25	0.71 mm	0.05 ( 90 %C)	0.1%	99.9%
n°35	0.5 mm	0.13 ( 80 %C)	0.2%	99.8%
n°60	0.25 mm	6.61 ( 0 %C)	6.8%	93.2%
n°80	0.18 mm	23.25 ( 0 %C)	30.0%	70.0%
n°120	0.125 mm	60.04 ( 0 %C)	90.1%	9.9%
n°230	0.063 mm	9.86 ( 0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.05 ( 0 %C)	100.0%	0.0%

% FINOS	0.05%
% CASCAJO	0.15%
<b>D16</b>	0.222
<b>D50</b>	0.162
<b>D84</b>	0.131
<b>Media</b>	0.171
<b>σ</b>	0.046
<b>M. distribuc.</b>	0.176

### MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm