


ANALISIS GRANULOMETRICOS

MUESTRA: D764

UTM X: 651013


UTM Y: 3220741



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	2.19 (100 %C)	2.2%	97.8%
n°18	1 mm	3.54 (100 %C)	5.7%	94.3%
n°25	0.71 mm	2.95 (100 %C)	8.7%	91.3%
n°35	0.5 mm	4.04 (0 %C)	12.7%	87.3%
n°60	0.25 mm	26.18 (0 %C)	38.9%	61.1%
n°80	0.18 mm	15.74 (0 %C)	54.6%	45.4%
n°120	0.125 mm	28.03 (0 %C)	82.7%	17.3%
n°230	0.063 mm	15.99 (0 %C)	98.7%	1.3%
Finos	<0,063 mm	1.35 (0 %C)	100.0%	0.0%

% FINOS	1.47%
% CASCAJO	8.67%
D16	0.399
D50	0.181
D84	0.114
Media	0.232
σ	0.142
M. distribuc.	0.257

MODA ARENAS FINAS



	<i>D50</i>	
<i>Cantos y gravas</i>	>2	mm
<i>Arena muy gruesa</i>	1-2	mm
<i>Arena gruesa</i>	0,5-1	mm
<i>Arena media</i>	0,25-0,5	mm
<i>Arena fina</i>	0,125-0,25	mm
<i>Arena muy fina</i>	0,062-0,125	mm
<i>Fango</i>	<0,062	mm