


## ANALISIS GRANULOMETRICOS

MUESTRA: IP132


UTM X: 622892

UTM Y: 3198032



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 ( 0 %C)	0.0%	100.0%
n°10	2 mm	0.00 ( 0 %C)	0.0%	100.0%
n°18	1 mm	0.07 ( 0 %C)	0.1%	99.9%
n°25	0.71 mm	0.64 ( 0 %C)	0.7%	99.3%
n°35	0.5 mm	4.11 ( 0 %C)	4.8%	95.2%
n°60	0.25 mm	55.71 ( 0 %C)	60.5%	39.5%
n°80	0.18 mm	31.04 ( 0 %C)	91.6%	8.4%
n°120	0.125 mm	7.58 ( 0 %C)	99.1%	0.9%
n°230	0.063 mm	0.85 ( 0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.00 ( 0 %C)	100.0%	0.0%

% FINOS	0.00%
% CASCAJO	0.00%
<b>D16</b>	0.450
<b>D50</b>	0.297
<b>D84</b>	0.197
<b>Media</b>	0.315
<b>σ</b>	0.126
<b>M. distribuc.</b>	0.323



### MODA ARENAS MEDIAS

	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm