


ANALISIS GRANULOMETRICOS

MUESTRA: D752

UTM X: 648313


UTM Y: 3242350



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.84 (100 %C)	0.8%	99.2%
n°10	2 mm	3.05 (100 %C)	3.9%	96.1%
n°18	1 mm	11.99 (100 %C)	15.9%	84.1%
n°25	0.71 mm	9.25 (100 %C)	25.1%	74.9%
n°35	0.5 mm	11.67 (0 %C)	36.8%	63.2%
n°60	0.25 mm	26.71 (0 %C)	63.5%	36.5%
n°80	0.18 mm	8.94 (0 %C)	72.5%	27.5%
n°120	0.125 mm	15.46 (0 %C)	87.9%	12.1%
n°230	0.063 mm	11.25 (0 %C)	99.2%	0.8%
Finos	<0,063 mm	0.84 (0 %C)	100.0%	0.0%

% FINOS	1.12%
% CASCAJO	25.13%
D16	0.497
D50	0.259
D84	0.124
Media	0.293
σ	0.186
M. distribuc.	0.311

MODA ARENAS MEDIAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm