


ANALISIS GRANULOMETRICOS

MUESTRA: IP133


UTM X: 622937

UTM Y: 3198122



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.00 (0 %C)	0.0%	100.0%
n°18	1 mm	0.06 (0 %C)	0.1%	99.9%
n°25	0.71 mm	0.69 (0 %C)	0.8%	99.2%
n°35	0.5 mm	5.02 (0 %C)	5.8%	94.2%
n°60	0.25 mm	65.95 (0 %C)	71.7%	28.3%
n°80	0.18 mm	23.14 (0 %C)	94.9%	5.1%
n°120	0.125 mm	5.02 (0 %C)	99.9%	0.1%
n°230	0.063 mm	0.06 (0 %C)	99.9%	0.1%
Finos	<0,063 mm	0.06 (0 %C)	100.0%	0.0%

% FINOS	0.06%
% CASCAJO	0.00%
D16	0.461
D50	0.332
D84	0.213
Media	0.335
σ	0.124
M. distribuc.	0.337



MODA ARENAS MEDIAS

	<i>D50</i>	
<i>Cantos y gravas</i>	>2	mm
<i>Arena muy gruesa</i>	1-2	mm
<i>Arena gruesa</i>	0,5-1	mm
<i>Arena media</i>	0,25-0,5	mm
<i>Arena fina</i>	0,125-0,25	mm
<i>Arena muy fina</i>	0,062-0,125	mm
<i>Fango</i>	<0,062	mm

