


ANALISIS GRANULOMETRICOS

MUESTRA: D731


UTM X: 646797

UTM Y: 3251421



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	6.92 (100 %C)	6.9%	93.1%
n°10	2 mm	18.79 (100 %C)	25.7%	74.3%
n°18	1 mm	26.95 (100 %C)	52.7%	47.3%
n°25	0.71 mm	8.28 (100 %C)	60.9%	39.1%
n°35	0.5 mm	7.42 (0 %C)	68.4%	31.6%
n°60	0.25 mm	14.83 (0 %C)	83.2%	16.8%
n°80	0.18 mm	6.67 (0 %C)	89.9%	10.1%
n°120	0.125 mm	6.80 (0 %C)	96.7%	3.3%
n°230	0.063 mm	2.60 (0 %C)	99.3%	0.7%
Finos	<0,063 mm	0.74 (0 %C)	100.0%	0.0%

% FINOS	1.90%
% CASCAJO	60.94%
D16	0.533
D50	0.296
D84	0.149
Media	0.326
σ	0.192
M. distribuc.	0.341



MODA ARENAS MEDIAS

	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm