


ANALISIS GRANULOMETRICOS

MUESTRA: D754

UTM X: 654594


UTM Y: 3225891



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.10 (100 %C)	0.1%	99.9%
n°10	2 mm	5.95 (100 %C)	6.1%	93.9%
n°18	1 mm	20.62 (100 %C)	26.7%	73.3%
n°25	0.71 mm	11.59 (100 %C)	38.3%	61.7%
n°35	0.5 mm	8.51 (0 %C)	46.8%	53.2%
n°60	0.25 mm	29.54 (0 %C)	76.3%	23.7%
n°80	0.18 mm	14.05 (0 %C)	90.4%	9.6%
n°120	0.125 mm	8.21 (0 %C)	98.6%	1.4%
n°230	0.063 mm	1.13 (0 %C)	99.7%	0.3%
Finos	<0,063 mm	0.31 (0 %C)	100.0%	0.0%

% FINOS	0.50%
% CASCAJO	38.26%
D16	0.488
D50	0.311
D84	0.181
Media	0.327
σ	0.154
M. distribuc.	0.335

MODA ARENAS MEDIAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm