


## ANALISIS GRANULOMETRICOS

MUESTRA: D760

UTM X: 651516


UTM Y: 3222830



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 ( 0 %C)	0.0%	100.0%
n°10	2 mm	1.20 ( 100 %C)	1.2%	98.8%
n°18	1 mm	10.60 ( 100 %C)	11.8%	88.2%
n°25	0.71 mm	9.64 ( 100 %C)	21.4%	78.6%
n°35	0.5 mm	15.18 ( 0 %C)	36.6%	63.4%
n°60	0.25 mm	28.80 ( 0 %C)	65.4%	34.6%
n°80	0.18 mm	10.84 ( 0 %C)	76.3%	23.7%
n°120	0.125 mm	16.02 ( 0 %C)	92.3%	7.7%
n°230	0.063 mm	7.11 ( 0 %C)	99.4%	0.6%
Finos	<0,063 mm	0.60 ( 0 %C)	100.0%	0.0%

% FINOS	0.77%
% CASCAJO	21.45%
<b>D16</b>	0.536
<b>D50</b>	0.291
<b>D84</b>	0.142
<b>Media</b>	0.323
<b>σ</b>	0.197
<b>M. distribuc.</b>	0.339

## MODA ARENAS MEDIAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm