


ANALISIS GRANULOMETRICOS

MUESTRA: D763

UTM X: 651482


UTM Y: 3221846



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.20 (100 %C)	0.2%	99.8%
n°10	2 mm	3.13 (100 %C)	3.3%	96.7%
n°18	1 mm	10.81 (100 %C)	14.1%	85.9%
n°25	0.71 mm	9.09 (100 %C)	23.2%	76.8%
n°35	0.5 mm	9.39 (0 %C)	32.6%	67.4%
n°60	0.25 mm	26.87 (0 %C)	59.5%	40.5%
n°80	0.18 mm	10.40 (0 %C)	69.9%	30.1%
n°120	0.125 mm	17.07 (0 %C)	87.0%	13.0%
n°230	0.063 mm	12.22 (0 %C)	99.2%	0.8%
Finos	<0,063 mm	0.81 (0 %C)	100.0%	0.0%

% FINOS	1.05%
% CASCAJO	23.23%
D16	0.473
D50	0.236
D84	0.121
Media	0.277
σ	0.176
M. distribuc.	0.297

MODA
ARENAS FINAS



	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm