


ANALISIS GRANULOMETRICOS

MUESTRA: IP63


UTM X: 647406

UTM Y: 3208887



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.07 (50 %C)	0.1%	99.9%
n°18	1 mm	0.21 (50 %C)	0.3%	99.7%
n°25	0.71 mm	0.70 (25 %C)	1.0%	99.0%
n°35	0.5 mm	5.21 (0 %C)	6.2%	93.8%
n°60	0.25 mm	73.87 (0 %C)	80.1%	19.9%
n°80	0.18 mm	17.11 (0 %C)	97.2%	2.8%
n°120	0.125 mm	2.75 (0 %C)	99.9%	0.1%
n°230	0.063 mm	0.07 (0 %C)	100.0%	0.0%
Finos	<0,063 mm	0.00 (0 %C)	100.0%	0.0%

% FINOS	0.00%
% CASCAJO	0.32%
D16	0.466
D50	0.351
D84	0.234
Media	0.350
σ	0.116
M. distribuc.	0.350



MODA ARENAS MEDIAS

	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm