


ANALISIS GRANULOMETRICOS

MUESTRA: D737

UTM X: 646466


UTM Y: 3249898



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	7.44 (100 %C)	7.4%	92.6%
n°10	2 mm	7.09 (100 %C)	14.5%	85.5%
n°18	1 mm	12.29 (100 %C)	26.8%	73.2%
n°25	0.71 mm	5.65 (100 %C)	32.5%	67.5%
n°35	0.5 mm	4.48 (0 %C)	37.0%	63.0%
n°60	0.25 mm	21.88 (0 %C)	58.8%	41.2%
n°80	0.18 mm	17.85 (0 %C)	76.7%	23.3%
n°120	0.125 mm	17.76 (0 %C)	94.4%	5.6%
n°230	0.063 mm	4.75 (0 %C)	99.2%	0.8%
Finos	<0,063 mm	0.81 (0 %C)	100.0%	0.0%

% FINOS	1.20%
% CASCAJO	32.47%
D16	0.428
D50	0.221
D84	0.141
Media	0.263
σ	0.143
M. distribuc.	0.285

MODA ARENAS FINAS



	<u>D50</u>	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm