


ANALISIS GRANULOMETRICOS

MUESTRA: D733

UTM X: 647239


UTM Y: 3251767



TAMICES ASTM		Peso retenido (gramos)	% Retenido	% Que pasa
n°4	4.75 mm	0.00 (0 %C)	0.0%	100.0%
n°10	2 mm	0.25 (100 %C)	0.3%	99.7%
n°18	1 mm	0.59 (100 %C)	0.8%	99.2%
n°25	0.71 mm	2.46 (100 %C)	3.3%	96.7%
n°35	0.5 mm	10.59 (0 %C)	13.9%	86.1%
n°60	0.25 mm	60.85 (0 %C)	74.7%	25.3%
n°80	0.18 mm	15.68 (0 %C)	90.4%	9.6%
n°120	0.125 mm	8.39 (0 %C)	98.8%	1.2%
n°230	0.063 mm	1.02 (0 %C)	99.8%	0.2%
Finos	<0,063 mm	0.17 (0 %C)	100.0%	0.0%

% FINOS	0.18%
% CASCAJO	3.31%
D16	0.480
D50	0.345
D84	0.206
Media	0.344
σ	0.137
M. distribuc.	0.343

MODA ARENAS MEDIAS



	D50	
Cantos y gravas	>2	mm
Arena muy gruesa	1-2	mm
Arena gruesa	0,5-1	mm
Arena media	0,25-0,5	mm
Arena fina	0,125-0,25	mm
Arena muy fina	0,062-0,125	mm
Fango	<0,062	mm